

Appendix I







ISO 15223-1:2007(E)

5 Symbols

When appropriate, information essential for proper use shall be indicated on the medical device, on its package or in the associated documentation by using the corresponding symbols given in Table 1.

Examples can be found in Annex A.

Table 1 — Symbols to convey information essential for proper use

No.	Symbol	Title	ISO 7000 or IEC 60417 registration number
5.1		Biological risks	ISO 7000-0659 (DB 2004-01)
5.2		Do not re-use	ISO 7000-1051 (DB 2004-01)
5.3		Consult instructions for use NOTE This symbol advises the reader to consult the operating instructions for information needed for the proper use of the device. See also symbol 5.4.	ISO 7000-1641 (DB 2004-01)
5.4		Caution, consult accompanying documents NOTE 1 This symbol advises the reader to consult the accompanying documents for important safety-related information such as warnings and precautions that cannot, for a variety of reasons, be presented on the device itself. See also symbol 5.3. NOTE 2 The symbol A or B in ISO 7000-0434 ("Caution") may also be used.	ISO 7000-0434 (DB 2004-01)
5.5		Fragile, handle with care	ISO 7000-0621 (DB 2004-01)
5.6		Keep away from sunlight NOTE The symbol may also mean "Keep away from heat" as referenced in ISO 7000:1989.	ISO 7000-0624 (DB 2004-01)

Source: ISO 15223-1 Medical devices – symbols to be used with medical device labels, labelling and information to be supplied – Part 1: general requirements